




Suggested Cataract Set



This set of instruments is ideal for the phacoemulsification cataract surgeon who prefers to perform surgical incisions using diamond knives and needs the facility to insert intra-ocular lens implants using a folding technique.


9-556 Kratz Barraquer Speculum

- 14.5mm open blades
- Angled to rest temporarily
- Single piece construction



9-588-3 Lieberman Adjustable Speculum


- 15mm open blades
- Angled to rest temporarily
- Adjustable with thumb screw



Blades when closed




2-100 Pierse Notched Forceps



- Pierse 0.25mm notched, 6mm tying platforms
- Straight shafts
- Standard handle, length 85mm




2-110N Straight Suturing forceps (0.12mm)



- 0.12mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- Standard handle, length 84mm




2-132 DK Troutman Barraquer Colibri Forceps



- 0.12mm, 1 x 2 teeth, tip length 2mm
- 6mm tying platforms
- Colibri style shafts
- Standard handle, length 84mm




2-501N Curved Tying Forceps



- 6mm tying platforms
- Curved shafts
- Standard handle, length 83mm



4-440 DK Angled Retractable Diamond Knife, 3mm Spear




W=3 T=0.17 L=4

- 3 spear diamond blade
- 45° angled retractable handle
- Length 135mm



4-590 DK Retractable Diamond Knife, 45° Double Edge




W=1 T=0.17 L=3.5

- 1mm wide 45° double edge diamond blade
- Retractable long handle
- Length 123.5mm



2-2-716G-9R Calladine-Inamura Round Handle Capsulorhexis Forceps, Corneal Incision





- Pointed serrated interlocking tips
- Curved shaft, tip to pivot point 10mm
- Marks on shaft at 2.5mm and 5mm denote desired size of capsulorhexis

- Cross action tips, 1.5mm width at pivot box
- Tips angled 45° to handle
- Round handle, overall length 120.5mm




6-085 Barrett Duo Nucleus Rotator / Manipulator / Splitter



Nucleus splitter


- Straight, sharp inner / outer sides of tip, tip length 1.25mm
- Cutting Edge 90° to Axis
- 45° angled shaft, tip to angle length 10mm



Rotator / manipulator

- 0.65mm mushroom tip
- 45° angled shaft, tip to angle length 10mm
- Barrett balanced set handle, length 123mm

8-652 / 8-652S Irrigation Handpiece's



- 23 gauge
- 0.65mm tube diameter
- Two 0.4mm irrigation ports
- Curved shaft, tube length 15mm
- Round handle, length 107mm (8-652)
- Round handle, length 80mm (8-652S)




8-652




8-652S

8-657 / 8-657S Aspiration Handpiece's



- 23 gauge
- 0.65mm tube diameter
- 0.3mm aspiration port
- Curved shaft, tube length 11mm
- Round handle, length 103mm (8-657)
- Round handle, length 77mm (8-657S)



8-657



8-657S


6-510 Mackool Capsule Polisher




- 0.3mm diameter olive shaped textured tip
- 45° angled shaft, tip to angle length 11mm
- Round handle, length 126mm




2-700 Shepard IOL Forceps



- Slightly curved tips
- Standard cross action handle, length 108mm



2-774-1 DK IOL Folding Forceps




- Paddle style jaws
- Stops in handle design prevent "over folding"
- Standard long handle, length 110mm



Conductive for upward folding. Highly polished inner jaw surface protects from scratching. Designed specifically for AcrySof approved by Alcon.

DK7740 DK IOL Insertion Forceps




- Biconvex jaw design
- 40° angled shafts, tip to angle length 7.5mm
- Standard long handle, length 110mm

6-450 Lester IOL Manipulator




- 0.4mm rounded tip / rotator, shaft to rotator length 0.8mm
- Push, pull style
- Angled shaft, tip to angle length 9mm
- Round handle, length 118mm

3-205 DK Needle Holder (curved)



- 7mm delicate curved jaws
- Round handle, tag spring
- Length 107mm

3-206 DK Needle Holder (curved with lock)



- 7mm delicate curved jaws
- Round handle with lock, tag spring
- Length 107mm


Sterilising Tray (suitable for 10 instruments - T7010 or 20 Instruments - T7020)

T7010

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat

T7020

- External dimensions: 39.5cm x 26.8cm x 3.2cm
- Internal dimensions: 38cm x 25.3cm x 1.5cm (without silicone mat)
- Suitable for 20 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© November 2014 Duckworth & Kent
07.07.2016

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeon's names or by the category of the instrument.

at the Leading Edge