



Suggested Comprehensive Lid Set



A comprehensive set of instruments allowing the surgeon to perform simple epilation of lashes and curettage of chalazia, to ectropion repairs and ptosis procedures.

2-640 DK Desmarres Chalazion Forceps / Clamp

- 11.5mm x 17.5mm inside oval circular working diameter upper plate
- 14mm x 20mm diameter lower plate
- Polished finish lower plate
- Standard handle, locking thumb screw, length 88mm




6-641 DK Meyerhoefer Chalazion Curette

- Straight shaft
- Round Handle,
- length 125.5mm




6-641




1.0mm

6-641-1




1.5mm

6-641-2




2.0mm

6-641-3




2.5mm

6-641-4




3.0mm

2-635 Kudo Meibomian Gland Pressing Forceps




Enlarged view on arrow



- 60° angled shafts, tip to angle length 11.5mm
- Internal faces of tip are cross hatched
- Standard cross action handle, length 112mm


2-2195 Cilia Forceps



- Tapered, smooth tips
- 30° angled shafts, tip to angle length 7mm
- Short ridged handle, length 85mm
- Single piece construction




2-401NR8 *Straight 0.65mm Notched Round Handle Forceps*





- Pierce 0.65mm notched, 6mm tying platforms
- Straight shafts
- 8mm diameter round handle, length 115mm

2-116N *Straight Suturing Forceps*




- 0.3mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- New long flat handle, length 115mm

1-500 *Westcott Style Stitch Scissors (Straight)*




- Semi-sharp tips, straight blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 115mm

1-512 *Westcott Style Tenotomy Scissors (Right Handed)*




- Semi-sharp tips, straight blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 115mm

6-805 *Bulldog Clip*




- 8.5mm serrated cross action jaws
- Length 46mm

6-628 Desmarres Lid Retractor, size 0



- Blade width 11mm
- Flat handle, length 134mm

6-629 Rabkin Lid Retractor



- Blade width 15.5mm
- Flat handle, length 134mm
- Non-reflective finish, textured surface top face

Lower lid retractor designed to function as a broad backstop in carbon dioxide laser transconjunctival blepharoplasty.

6-664 Otaka Lid Plate

- 23.5mm wide
- 20.4mm wide
- Length 55mm




6-665 Rabkin Laser Blepharoplasty Lid Plate

- 20mm wide
- 24mm wide
- Inside highly polished
- Partial matte finish
- Length 111mm
- Non-reflective finish

For use in carbon dioxide laser upper lid and transconjunctival blepharoplasty as well as to retract eyelashes in laser resurfacing of skin.




6-625 Jameson Muscle Hook




- 2mm bulbous tip, tip to curve length 8mm
- Round handle, length 135mm

6-620 DK Paediatric Muscle Hook



- 0.65mm diameter tip, tip to curve length 4.5mm
 - Round handle, length 125 vmm
- Used in inferior oblique surgery


9-650 DK Castroviejo Style Marking Caliper



- 0.25mm x 0.25mm delicate marking tips
- Marks 0mm to 15mm in 0.25 increments
- Marks from centre of tips
- Adjustable thumb screw
- Standard caliper handle




3-202 DK Needle Holder (curved with lock)



- 9mm medium curved jaws
- Round handle with lock, tag spring
- Length 109mm



3-203 DK Needle Holder (curved)




- 9mm medium curved jaws
- Round handle, tag spring
- Length 109mm



T7010 Sterillising Tray (suitable for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© August 2015 Duckworth & Kent
04.07.2016

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeon's names or by the category of the instrument.

at the Leading Edge