



## Suggested Squint Set




A comprehensive set of instruments that provide the squint surgeon with the necessary specula, hooks, forceps, scissors and needle holders.

### 9-599 *Shepard Reversible Speculum*

- 14mm solid blades
- Angled to rest temporally or nasally
- Adjustable with thumb screw




### 2-501-2 *Chihara Curved Conjunctival Forceps*




- 6mm tying platforms
- Serrated tips
- Curved shafts
- Standard handle, length 88mm




### 2-510-1N *Moorfields Utility Forceps*




- Delicate serrated tips, 1mm wide
- Straight shafts
- New long flat handle, length 115mm




### 1-501 *Westcott Style Stitch Scissors (curved)*




- Semi-sharp tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm




### 1-510 Westcott Style Tenotomy Scissors (curved)




- Blunt tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm




### 6-628 Desmarres Lid Retractor (size 0)



- Blade width 11mm
- Flat handle, length 134mm




### 6-620 DK Paediatric Muscle Hook



- 0.65mm diameter tip, tip to curve length 4.5mm
- Flat handle, length 125mm



### 6-625 Jameson Muscle Hook




- 2mm bulbous tip, tip to curve length 8mm
- Flat handle, length 135mm



### 6-626 Assaf Resection Muscle Hook


Also available with hook sizes of 12mm (ref: 6-626-1) and 14mm (ref: 6-626-2)



- 90° double 10mm adjustable tips
- Caliper measures 2mm - 15mm
- Flat handle, length 135mm




### 9-650 DK Castroviejo Style Marking Caliper




- 0.25mm x 0.25mm delicate marking tips
- Marks 0mm to 15mm in 0.25 increments
- Marks from centre of tips
- Adjustable thumb screw
- Standard caliper handle




### 1-500 Westcott Style Stitch Scissors (straight)



- Semi-sharp tips, straight blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm




### 1-500B DK Westcott Style Stitch Scissors (straight)



- Blunt tips, straight blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 115mm




### 2-114N Straight Toothed Forceps (0.12mm) - Flat Handle



- 0.12mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- New long flat handle, length 115mm




### 2-116N Straight Toothed Forceps (0.3mm)



- 0.3mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- New long flat handle, length 115mm




### 2-504N Straight Tying Forceps




- 6mm tying platforms
- Straight shafts
- New long flat handle, length 115mm

### 3-202 DK Needle Holder (curved with lock)



- 9.0mm medium curved jaws
- Round handle with lock, tag spring
- Length 109mm


### 3-203 DK Needle Holder (curved)



- 9.0mm medium curved jaws
- Round handle, tag spring
- Length 109mm

### T7020 20 Instrument Sterilising Tray

- External dimensions: 39.5cm x 26.8cm x 3.2cm
- Internal dimensions: 38cm x 25.3cm x 1.5cm (without silicone mat)
- Suitable for 20 instruments
- One silicone mat



[www.duckworth-and-kent.com](http://www.duckworth-and-kent.com)

7 Marquis Business Centre  
Royston Road, Baldock  
Herts SG7 6XL England

Tel: +44 (0)1462 893254  
Fax: +44 (0)1462 896288  
Email: [info@duckworth-and-kent.co.uk](mailto:info@duckworth-and-kent.co.uk)



© March 2011 Duckworth & Kent  
06.12.2016

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

*at the Leading Edge*