



Suggested Basic Lid Set




This set provides the surgeon with an array of instruments that will allow for a variety of eyelid surgical procedures including the correction of ectropia and entropia, as well as ptosis repairs and eyelid trauma. The set also includes a curette and clamp for removing chalazia.

6-664 *Otaka Lid Plate*

- 23.5mm wide
- 20.4mm wide
- Length 55mm




2-401NR8 *Straight 0.65mm Notched Round Handle Forceps*





- Pierce 0.65mm notched, 6mm tying platforms
- Straight shafts
- 8mm diameter round handle, length 115mm

2-116N *Straight Toothed Forceps (0.3mm)*




- 0.3mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- New long flat handle, length 115mm

1-500 *Westcott Style Stitch Scissors (straight)*




- Semi-sharp tips, straight blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm


1-512 Westcott Style Tenotomy Scissors (Right Handed)



- Sharp tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm




6-805 Bulldog Clip



- 8.5mm serrated cross action jaws
- Length 46mm




6-628 Desmarres Lid Retractor (size 0)



- Blade width 11mm
- Flat handle, length 134mm




3-202 DK Needle Holder (curved with lock)



- 9mm medium curved jaws
- Round handle with lock, tag spring
- Length 109mm



3-203 DK Needle Holder (curved)



- 9mm medium curved jaws
- Round handle tag spring
- Length 109mm




2-640 DK Desmarres Chalazion Forceps / Clamp

- 11.5mm x 17.5mm inside oval circular working diameter upper plate
- 14.0mm x 20mm diameter lower plate
- Polished finish lower plate
- Standard handle, locking thumb screw, length 88mm




6-641 DK Meyerhoefer Chalazion Curette, 1mm cup



- 1mm cup
- Straight shaft
- Round Handle, length 125.5mm

T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



Visit : www.duckworth-and-kent.com to view
over 1000 ophthalmic titanium instruments

The screenshot shows the Duckworth & Kent Ltd. website. The header includes the company logo and name, with navigation links for 'about', 'services', 'business partners', 'careers', and 'contact'. Below the header is a main navigation bar with 'Home', 'Products', 'Events', and 'Literature/Videos'. A search bar is located on the right. The main content area features a Facebook social media link, a description of the company as a world leader in titanium ophthalmic surgical instrumentation, and a 'FEATURE PRODUCTS' section. This section includes images and names of products such as 'Toric Markers', 'Femtosecond Laser Instruments', 'Capsulorhexis Range', 'Femtosecond Laser', 'Toric Marker', 'N Forceps', and 'AMO Injector'. A 'Connect with Duckworth & Kent' section displays icons for Facebook, LinkedIn, YouTube, and Google+. A QR code is positioned to the left of the website screenshot. To the right, there is a globe icon and the slogan 'at the Leading Edge'. The footer of the website screenshot provides contact information: Tel: +44 (0)1462 893254, Fax: +44 (0)1462 896288, and email: info@duckworth-and-kent.co.uk. It also includes a secondary navigation bar and copyright information for 2014.

As well as product information for over 1000 instruments, you will also find detailed literature and informative videos showing some of our instruments in action.

Connect with Duckworth & Kent



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© November 2014 Duckworth & Kent
04.07.2016

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeons' names or by the category of the instrument.

at the Leading Edge