



Suggested Vitreoretinal Set




6-676 DK Squeeze Handle for VR Instrument Heads

- Suitable for all Duckworth & Kent Vitreoretinal Instrument Heads
- Overall Length 92 mm (without head)
- Squeeze action handle



Scissors



- Sharp tipped straight micro scissors with 2.2mm long blade
- Straight 23 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

1-841N



23 Gauge



- Blunt tipped straight micro scissors with 1.9mm long blade
- Straight 23 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

1-841BN



23 Gauge



- Sharp tipped curved micro scissors with 2.2mm long blade
- Straight 23 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

1-842N



23 Gauge



- Sharp tipped straight micro scissors with 2.5mm long blade
- Straight 20 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

1-841-1N



20 Gauge

Forceps



End Gripping Forceps
ILM Peeling

23G

- Fine end gripping tips, tip length 3.2mm
- Straight 23 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-877N



23 Gauge

Delicate Serrated
Forceps

23G

- Delicate serrated jaws with blunt tip, tip length 3.0mm
- Straight 23 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-878-1N



23 Gauge

Asymmetrical Forceps

NEW tip design 23G

- Asymmetrical tip design with fine platform
- Straight 23 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-879N



23 Gauge

Rassam Membrane
Picking Forceps

20G

- 2.0mm jaws, 0.7mm platforms, 55° angled pointed picking tip
- Straight 20 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-885N



20 Gauge

Smooth Rounded Tip
Forceps

20G

- 2.0mm smooth rounded tips
- Straight 20 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-868N



20 Gauge

End Gripping Forceps
ILM Peeling

20G

- Fine end gripping tips, tip length 3.7mm
- Straight 20 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-871N



20 Gauge

Delicate Serrated
Forceps

20G

- Delicate serrated jaws with blunt tip, tip length 3.0mm
- Straight 20 gauge tube, length 32.5mm
- Requires DK Squeeze Handle, ref 6-676
- Detachable tip heads for cleaning purposes
- Once attached to round squeeze handle, length is 143mm

2-878N



20 Gauge

Fixed Heads

20 and 23 Gauge Forceps & Scissors

The fixed head range of Vitreoretinal forceps and scissors are a one piece devices featuring the same tip and handle design as the interchangeable head.

All fixed head vitreoretinal forceps and scissors are supplied with a protective guard which can be sterilised.



Scissors



Scissors Variations




Straight Scissors

23G

- Sharp tipped straight micro scissors with 2.2mm long blade
- Straight 23 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 138mm

1-841



23 Gauge




Straight Blunt Scissors

23G

- Blunt tipped straight micro scissors with 1.75mm long blade
- Straight 23 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 138mm

1-841B



23 Gauge




Curved Scissors

NEW tip design 23G

- Sharp tipped curved micro scissors with 2.5mm long blade
- Straight 23 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 138mm

1-842



23 Gauge




Straight Scissors

20G

- Sharp tipped straight micro scissors with 2.5mm long blade
- Straight 20 gauge tube, length 32.5mm
- Squeeze action handle closes scissors
- Round squeeze handle, length 138mm

1-841-1



20 Gauge

Duckworth & Kent - Titanium

Forceps



Forceps Variations




End Gripping Forceps
ILM Peeling

23G

- Fine end gripping tips, tip length 3mm
- Straight 23 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 138mm

2-877



23 Gauge




Delicate Serrated
Forceps

23G

- Delicate serrated jaws with blunt tip, tip length 3mm
- Straight 23 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 139mm

2-878-1



23 Gauge




Asymmetrical Forceps

NEW tip design 23G

- Asymmetrical tip design with fine platforms
- Straight 23 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 138mm

2-879



23 Gauge




Rassam Membrane
Picking Forceps

20G

- 2.0mm jaws, 0.7mm platforms, 55° angled pointed picking tips
- Straight 20 gauge tube, tube length 32.5mm
- Squeeze action activates both jaws
- Round squeeze handle, length 138mm

2-885



20 Gauge




Smooth Rounded Tip
Forceps

20G

- 2.0mm smooth rounded tips
- Straight 20 gauge tube, tube length 32.5mm
- Squeeze action activates both jaws
- Round squeeze handle, length 138mm

2-868



20 Gauge




End Gripping Forceps
ILM Peeling

20G

- Fine end gripping tips, tip length 3.7mm
- Straight 20 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 141mm

2-871



20 Gauge




Delicate Serrated
Forceps

20G


- Delicate serrated jaws with blunt tip, tip length 3mm
- Straight 20 gauge tube, tube length 32.5mm
- Squeeze handle activates both jaws
- Round squeeze handle, overall length 138mm

2-878



20 Gauge


6-800 Barrett LeClip Utility Clamp



- 14mm serrated cross action jaws
- Length 82mm
- Distinctive identification labelling




6-805 Bulldog Clip



- 8.5mm serrated cross action jaws
- Length 46mm




9-650 DK Castroviejo Style Marking Caliper



- 0.25mm x 0.25mm delicate marking tips
- Marks 0.0mm to 15mm in 0.25 increments
- Marks from centre of tips
- Adjustable thumb screw
- Standard caliper handle




2-160 Troutman Superior Rectus Forceps



- 0.5mm, 1 x 2 teeth
- 45° angled shafts, tip to angle length 10mm
- Standard long handle, length 111mm




2-118 DK Bonn Suturing Forceps (0.30mm)



- 0.3mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- Standard handle, length 88mm




2-132N Colibri Toothed Forceps (0.12mm)



- 0.12mm, 1 x 2 teeth, tip length 2mm
- 6mm tying platforms
- Colibri style shafts
- New flat handle, length 83mm




2-500N Straight Tying Forceps




- 6mm tying platforms
- Straight shafts
- New flat handle, length 84mm




1-118 DK Castroviejo Miniature Corneal Scissors (curved)




- Blunt tips, curved blades
- Universal cutting blades
- Cut length 6mm, tip to pivot length 11mm
- Standard handle, length 116mm




1-510 Westcott Style Tenotomy Scissors (curved)



- Blunt tips, curved blades
- Cut length 12mm, blade length 20mm
- Standard handle, length 116mm



4-600 Pearce Retractable Diamond Knife, 45° Double Edge



W=1.0T=0.17L=3.5

- Retractable long handle with second position 10mm shaft extension
- Length 123.5mm, extended length 133.5mm



6-912 *Giunchiglia Membrane Peeling Spatula, 23 Gauge*



- Ultra fine textured tip
- 0.60mm wide, blunt tip curved 45°



- Straight 23 gauge shaft, length 32.5mm
- Round handle, overall length 128mm

6-912-1 *Membrane Peeling Spatula, 23 Gauge*




- Fine textured tip
- 0.60mm wide, blunt tip curved 45°





- Straight 23 gauge shaft, length 32.5mm
- Round handle, overall length 128mm

9-556 *Kratz Barraquer Speculum*

- 14.5mm open blades
- Angled to rest temporarily
- Single piece construction




6-635-2 *DK Schocket Double Ended Scleral Depressor*



- 2.5mm round and 4.5mm cylinder tips
- Round handle, length 135mm




6-625 *Jameson Muscle Hook*



- 2mm bulbous tip, tip to curve length 8mm
- Flat handle, length 135mm




3-302 DK Barraquer Needle Holder (curved with lock)



- 9mm delicate curved jaws
- Round handle with lock, hinge spring
- Length 109mm



3-203 DK Needle Holder (curved)






- 9.0mm medium curved jaws
- Round handle, tag spring
- Length 109mm




6-675 Protective Cleaning Cover

- Used to protect all detachable Vitreo-Retinal tip heads
- Luer flushing adaptor



T7020 20 Instrument Sterilising Tray

- External dimensions: 39.5cm x 26.8cm x 3.2cm
- Internal dimensions: 38cm x 25.3cm x 1.5cm (without silicone mat)
- Suitable for 20 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© April 2011 Duckworth & Kent
30.09.16

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeon's names or by the category of the instrument.

at the Leading Edge