



Suggested Laser Blepharoplasty Set



A set of instruments to assist conventional and laser blepharoplasty surgical procedure's

6-667-4 Egi Rabkin Small Eye Shield

- 18mm wide x 20mm long
- Inside highly polished
- Supplied as one (1) each



6-667 Rabkin Eye Shield

- 20.5mm wide x 22mm long
- Inside highly polished
- Supplied as one (1) each



6-667-3 Egi Rabkin Medium Eye Shield

- 18.5mm wide x 21mm long
- Inside highly polished
- Supplied as one (1) each



6-667-2 Egi Rabkin Large Eye Shield

- 19mm wide x 22mm long
- Inside highly polished
- Supplied as one (1) each



6-665 Rabkin Laser Blepharoplasty Plate

- 20mm wide
- 24mm wide
- Inside highly polished
- Partial matte finish
- Length 111mm
- Non-reflective finish



6-664 Otaka Lid Plate

- 23.5mm wide
- 20.4mm wide
- Length 55mm



2-645 Rabkin Lid Clamp

- Clamp length 24mm
- Plate dimensions 31mm x 17mm
- Flat handle, length 103mm
- Non-reflective finish, textured surface upper top face




9-599 Shepard Reversible Speculum

- 14mm solid blades
- Angled to rest temporally or nasally





6-805 *Bulldog Clip*




- 8.5mm serrated cross action jaws
- Length 46mm

6-629 *Rabkin Lid Retractor*




- Blade width 15.5mm
- Flat handle, length 134mm
- Non-reflective finish, textured surface top face

2-401NR8 *Straight 0.65mm Notched Round Handle Forceps*




- Pierce 0.65mm notched, 6mm tying platforms
- Straight shafts
- 8mm diameter round handle, length 115mm

33-225 *DK Needle Holder (curved)*



- 12mm heavy curved jaws
- Round long handle, tag spring
- Length 139mm


2-108N Straight Suturing Forceps (0.50mm)



- 0.5mm, 1 x 2 teeth, 6mm tying platforms
- Straight shafts
- New long flat handle, length 116mm

T7010 Sterilising Tray (for 10 instruments)

- External dimensions: 26.4cm x 16.2cm x 2.6cm
- Internal dimensions: 25.4cm x 15.2cm x 1.8cm (without silicone mat)
- Suitable for 10 instruments
- One silicone mat



www.duckworth-and-kent.com

7 Marquis Business Centre
Royston Road, Baldock
Herts SG7 6XL England

Tel: +44 (0)1462 893254
Fax: +44 (0)1462 896288
Email: info@duckworth-and-kent.co.uk



© March 2011 Duckworth & Kent
04.07.2016

D&K® is a registered trademark. All other brand names are trademarks or registered trademarks of their respective owners. All schematic line drawings, photographs and copy in this leaflet are fully protected by copyright. No part of this leaflet may be reproduced in any form without prior written permission. We reserve the right to make changes at any time, without notice, in product specifications and availability. Descriptive, typographic, or photographic errors are subject to correction. Name(s) of instruments are often comprised of surgeon's name, combination of surgeon's names or by the category of the instrument.

at the Leading Edge